

## .ins file (shelx input)

User input → \* TITL 514\_RTJ\_I\_113\_1 in P222 #16  
 \* REM reset to P222 #16  
 \* CELL 0.72932 9.889 10.073 51.938 90 90 90  
 ZERR 4 0.00198 0.00201 0.01039 0 0 0  
 \* LATT -1  
 SYMM -X,-Y,+Z  
 SYMM +X,-Y,-Z  
 SYMM -X,+Y,-Z  
 SFAC C N Sm  
 UNIT 288 48 8  
 TREF  
 HKLF 4  
 END|

WEB: International tables  
 for crystallography

User input

SMX defaults

\* Found in XSCALE.HKL

CELL	wavelength, UC
ZERR	Z ( <u>estimated</u> # formula units per cell), estimated error in UC dimensions – might be able to just use all zeros here – SHELXL refinement updates this anyway
LATT	lattice type - negative if the structure is non-centrosymmetric
SYMM	symmetry operator: X, Y, Z is assumed and thus omitted.
SFAC	atom names in order of scattering factors to be employed by program
UNIT	# of each type of atom in UC in SFAC order
TREF	tells SHELX I want to use direct methods
HKLF	"n" reflection file reading instructions n = 4 is general purpose for SMX n = 3 for MX IR or anomalous